

Function:

- 1-Vacuum
- 2-Spray
- 3-Galvanic
- 4-High-Frequency

Instruction:

1. Vacuum & Spray
1- Connect small suction glass tube to the suction/spray hose, put the hose to the suction port of machine, press the small hole on the small glass tube, then machine is ready for suction operation.
- 2- Connect spray bottles to the suction/spray hose, put the hose to the spray port of machine, press the small hole on spray bottle then the machine is ready for spray like mist.
2. Galvanic
1- Assemble the accessories and connect black plug of Cathode Stick, and put the red plug of Bead or Scroll Bar, users use the Cathode Stick, while controller hold the Bead or Scroll Bar.
- 2-Turn the power switch ON, ready for operation.
- 3-Press the switch for Galvanic.
- 4-Adjust for the energy of Galvanic, the energy lamp of the Galvanic will change at the same time.
3. High-Frequency
1- Insert the Ozone electrode into the handle and keep touch well.
- 2- Connect the plug to the power outlet then turn on the machine.
- 3- Turn the button of energy lever to choose appropriate level according to skin conditions and type then begin to operate.
- 4- After over the operation, turn the button of energy lever to "O" mark and turn off power, pull out the electrode and clean well and dry storage.

Caution:

- 1- There are three types of glasses of which should be used separately according to any specific parts of the body.
After over the operation, used glass tubes should be cleaned
- 2- and sterilized and dry storage.
- 3- Take off all metal decorations (earrings, necklace etc.) before operation.
- 4- During the operation, not to stay too long time in a treated area to avoid scars.
- 5- After over operation, used electrodes should be cleaned and dry storage.

Parts:

- Red and black tube: 1 Electrotherapy stick: 1
Electrotherapy tube: 1 Spray the bottle: 1
The ducting leads: 1 Glass straw: 1